

Product Datasheet

Dog TNF ALPHA protein (orb435755)

Description	Recombinant of dog TNF ALPHA protein
Reactivity	Canine
Tested Applications	ELISA, WB
Target	TNF ALPHA
Preservatives	10% Trehalose. Phosphate buffered saline
Form/Appearance	Purified recombinant protein - lyophilized
Storage	-20°C only (ship +4°C)
Note	For research use only
Clonality	Recombinant
Purity	Rec. Protein
Expiration Date	6 months from date of receipt.