

Product Datasheet

Human PDGF protein (orb435767)

Description	Recombinant of human PDGF protein
Reactivity	Human
Tested Applications	ELISA, FA, WB
Target	PDGF
Preservatives	None present. 30 mM Acetic Acid
Form/Appearance	Purified recombinant protein - lyophilized
Concentration	0.1 mg/ml
Storage	-20°C only (ship +4°C)
Note	For research use only
Application notes	BB HOMODIMER
Clonality	Recombinant
Purity	Rec. Protein
Expiration Date	6 months from date of receipt.