

Product Datasheet

Human IL10 protein (orb435777)

Description	Recombinant of human IL10 protein
Reactivity	Human
Tested Applications	ELISA, FA, WB
Target	IL-10
Preservatives	None present. Sodium Phosphate
Form/Appearance	Purified recombinant IL-10 - lyophilized
Concentration	0.1mg/ml
Storage	-20°C only (ship +4°C)
Note	For research use only
Clonality	Recombinant
Purity	Rec. Protein
Expiration Date	6 months from date of receipt.