

Product Datasheet

Human IFN GAMMA protein (orb435780)

Description	Recombinant of human IFN GAMMA protein
Reactivity	Human
Tested Applications	ELISA, FA, WB
Target	IFN GAMMA
Preservatives	None present. Sodium Phosphate
Form/Appearance	Recombinant protein produced in E. coli - lyophilized
Concentration	1.0 mg/ml
Storage	-20°C only (ship +4°C)
Note	For research use only
Clonality	Recombinant
Purity	Rec. Protein
Expiration Date	6 months from date of receipt.