

Product Datasheet

Human IL1RA protein (orb435791)

Description	Recombinant of human IL1RA protein
Reactivity	Human
Tested Applications	ELISA, FA, WB
Target	IL-1RA
Preservatives	None present. TRIS buffered saline
Form/Appearance	Purified recombinant protein - lyophilized
Concentration	0.1mg/ml
Storage	-20°C only (ship +4°C)
Note	For research use only
Clonality	Recombinant
Purity	Rec. Protein
Expiration Date	6 months from date of receipt.