

Product Datasheet

Human G-CSF protein (orb435792)

Description	Recombinant of human G-CSF protein
Reactivity	Human
Tested Applications	FA
Target	G-CSF
Preservatives	None present. Lyophilised from 5mM Sodium Citrate pH 4.0
Form/Appearance	Purified recombinant protein - lyophilized
Concentration	0.1 mg/ml
Storage	-20°C only (ship +4°C)
Note	For research use only
Clonality	Recombinant
Purity	Rec. Protein
Expiration Date	6 months from date of receipt.