

Product Datasheet

Human SDF-1 ALPHA protein (orb435808)

Description	Recombinant of human SDF-1 ALPHA protein
Reactivity	Human
Tested Applications	ELISA, FA, WB
Target	SDF-1 ALPHA
Preservatives	None present
Form/Appearance	Purified recombinant protein - lyophilized
Concentration	0.1 mg/ml
Storage	-20°C only (ship +4°C)
Note	For research use only
Clonality	Recombinant
Purity	Rec. Protein
Expiration Date	6 months from date of receipt.