

Product Datasheet

Human ALPHA FETOPROTEIN protein (orb435820)

Description	Human ALPHA FETOPROTEIN protein
Reactivity	Human
Tested Applications	ELISA, IA
Target	ALPHA FETOPROTEIN
Preservatives	0.1% Sodium Azide (NaN ₃). TRIS buffered saline
Form/Appearance	Purified human alpha fetoprotein - liquid
Concentration	1.0 mg/ml
Storage	+4°C DO NOT FREEZE
Note	For research use only
Purity	Purified
Expiration Date	6 months from date of receipt.