

Product Datasheet

Human GALECTIN-1 protein (orb435834)

Description	Recombinant of human GALECTIN-1 protein
Reactivity	Human
Tested Applications	ELISA, FA, WB
Target	GALECTIN-1
Preservatives	None present. 10mM Sodium phosphate pH 7.5
Form/Appearance	Purified recombinant protein - lyophilized
Concentration	0.1 mg/ml
Storage	-20°C only (ship +4°C)
Note	For research use only
Clonality	Recombinant
Purity	Rec. Protein
Expiration Date	6 months from date of receipt.