

## Product Datasheet

### Human ADIPONECTIN protein (orb435838)

|                            |                                                |
|----------------------------|------------------------------------------------|
| <b>Description</b>         | Recombinant of human ADIPONECTIN protein       |
| <b>Reactivity</b>          | Human                                          |
| <b>Tested Applications</b> | ELISA, FA                                      |
| <b>Target</b>              | ADIPONECTIN                                    |
| <b>Preservatives</b>       | None present. 10mM TRIS, pH 8.575mM L-Arginine |
| <b>Form/Appearance</b>     | Purified recombinant protein - lyophilized     |
| <b>Concentration</b>       | 1.0 mg/ml                                      |
| <b>Storage</b>             | -20°C only (ship +4°C)                         |
| <b>Note</b>                | For research use only                          |
| <b>Clonality</b>           | Recombinant                                    |
| <b>Purity</b>              | Rec. Protein                                   |
| <b>Expiration Date</b>     | 6 months from date of receipt.                 |