

Product Datasheet

Human PAI-1 protein (orb435840)

Description	Recombinant of human PAI-1 protein
Reactivity	Human
Tested Applications	FA
Target	PAI-1
Preservatives	None present. 10mM Sodium Phosphate
Form/Appearance	Purified recombinant protein - lyophilized
Concentration	0.1mg/ml
Storage	-20°C only (ship +4°C)
Note	For research use only
Clonality	Recombinant
Purity	Rec. Protein
Expiration Date	6 months from date of receipt.