

Product Datasheet

Human DEFENSIN ALPHA 1 protein (orb435846)

Description	Recombinant of human DEFENSIN ALPHA 1 protein
Reactivity	Human
Tested Applications	ELISA, FA, WB
Target	DEFENSIN ALPHA 1
Preservatives	None present
Form/Appearance	Purified recombinant protein - lyophilized
Concentration	0.1mg/ml
Storage	-20°C only (ship +4°C)
Note	For research use only
Clonality	Recombinant
Purity	Rec. Protein
Expiration Date	6 months from date of receipt.