

Product Datasheet

Human BDNF protein (orb435850)

Description	Recombinant of human BDNF protein
Reactivity	Human
Tested Applications	FA
Target	BDNF
Preservatives	Purified recombinant protein - lyophilized
Form/Appearance	Purified recombinant protein - lyophilized
Concentration	0.1mg/ml
Storage	-20°C only (ship +4°C)
Note	For research use only
Clonality	Recombinant
Purity	Rec. Protein
Uniprot ID	P23560
Expiration Date	6 months from date of receipt.