

Product Datasheet

Human CD120b protein [Out of stock] (orb435851)

Description	Recombinant of human CD120b protein
Reactivity	Human
Tested Applications	ELISA, FA, WB
Target	CD120b
Preservatives	None present. Phosphate buffered saline.
Form/Appearance	Purified recombinant protein - lyophilized
Concentration	0.1mg/ml
Storage	-20°C after reconstitution
Note	For research use only
Clonality	Recombinant
Purity	Rec. Protein
Expiration Date	6 months from date of receipt.