

Product Datasheet

Human ANNEXIN II protein (orb435855)

Description	Recombinant of human ANNEXIN II protein
Reactivity	Human
Tested Applications	ELISA, WB
Target	ANNEXIN II
Preservatives	None present. Phosphate buffered saline.
Form/Appearance	Purified recombinant protein - liquid
Concentration	0.1 mg/ml
Storage	-70°C (ship frozen)
Note	For research use only
Clonality	Recombinant
Purity	Rec. Protein
Expiration Date	6 months from date of receipt.