

## Product Datasheet

### Human CD141 protein (orb435874)

|                            |                                            |
|----------------------------|--------------------------------------------|
| <b>Description</b>         | Recombinant of human CD141 protein         |
| <b>Reactivity</b>          | Human                                      |
| <b>Tested Applications</b> | FA                                         |
| <b>Target</b>              | CD141                                      |
| <b>Preservatives</b>       | None present. 10mM Sodium Phosphate        |
| <b>Form/Appearance</b>     | Purified recombinant protein - lyophilized |
| <b>Concentration</b>       | 0.1mg/ml                                   |
| <b>Storage</b>             | -20°C after reconstitution                 |
| <b>Note</b>                | For research use only                      |
| <b>Clonality</b>           | Recombinant                                |
| <b>Purity</b>              | Rec. Protein                               |
| <b>Expiration Date</b>     | 6 months from date of receipt.             |