

Product Datasheet

Human MMP-2 protein (orb435878)

Description	Recombinant of human MMP-2 protein
Reactivity	Human
Tested Applications	FA
Target	MMP-2
Preservatives	None present. Sodium phosphate Calcium chloride
Form/Appearance	Purified recombinant protein - lyophilized
Concentration	0.1mg/ml
Storage	-20°C only (ship +4°C)
Note	For research use only
Clonality	Recombinant
Purity	Rec. Protein
Expiration Date	6 months from date of receipt.