

Product Datasheet

Human C-REACTIVE PROTEIN protein (orb435884)

Description	Native human C-REACTIVE PROTEIN protein
Reactivity	Human
Tested Applications	ELISA
Target	C-REACTIVE PROTEIN
Preservatives	0.1% Sodium Azide (NaN ₃). Tris buffered saline, 5 mM CaCl ₂ .
Form/Appearance	Purified protein - liquid
Concentration	1.0 mg/ml
Storage	+4°C DO NOT FREEZE
Note	For research use only
Purity	Purified
Expiration Date	6 months from date of receipt.