

Product Datasheet

Human KAPPA LIGHT CHAIN protein (orb435888)

Description	Native human KAPPA LIGHT CHAIN protein
Reactivity	Human
Tested Applications	ELISA
Target	KAPPA LIGHT CHAIN
Preservatives	None present. Phosphate buffered saline
Form/Appearance	Purified human kappa light chain - liquid
Concentration	5.0 mg/ml
Storage	-20°C only (ship +4°C)
Note	For research use only
Purity	Preservative Free
Expiration Date	6 months from date of receipt.