

Product Datasheet

Human LAMBDA LIGHT CHAIN protein (orb435889)

Description	Native human LAMBDA LIGHT CHAIN protein
Reactivity	Human
Tested Applications	ELISA
Target	LAMBDA LIGHT CHAIN
Preservatives	0.1% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Purified human Ig lambda light chain - liquid
Concentration	5.0 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Purity	Purified
Expiration Date	6 months from date of receipt.