

## Product Datasheet

### Human LAMBDA LIGHT CHAIN protein (orb435890)

<b>Description</b>	Native human LAMBDA LIGHT CHAIN protein
<b>Reactivity</b>	Human
<b>Tested Applications</b>	ELISA
<b>Target</b>	LAMBDA LIGHT CHAIN
<b>Preservatives</b>	Purified human Ig lambda light chain - liquid
<b>Form/Appearance</b>	Purified human Ig lambda light chain - liquid
<b>Concentration</b>	5.0 mg/ml
<b>Storage</b>	-20°C only (ship +4°C)
<b>Note</b>	For research use only
<b>Purity</b>	Preservative Free
<b>Uniprot ID</b>	<b>POCG04</b>
<b>Expiration Date</b>	6 months from date of receipt.