

Product Datasheet

PLASMODIUM FALCIPARUM HRP2 protein (orb435895)

Description	Recombinant of PLASMODIUM FALCIPARUM HRP2 protein
Tested Applications	ELISA
Target	PLASMODIUM FALCIPARUM HRP2
Preservatives	0.1% Sodium Azide (NaN ₃). Sodium carbonate buffer, pH 9.6
Form/Appearance	Purified recombinant protein - liquid
Concentration	1.0 mg/ml
Storage	-20°C only (ship frozen)
Note	For research use only
Clonality	Recombinant
Purity	Rec. Protein
Expiration Date	6 months from date of receipt.