

Product Datasheet

Viral HEPATITIS B SURFACE ANTIGEN AD protein (orb435896)

Description	Recombinant of viral HEPATITIS B SURFACE ANTIGEN AD protein
Reactivity	Virus
Tested Applications	ELISA, WB
Target	HEPATITIS B SURFACE ANTIGEN AD
Preservatives	0.1% Sodium Azide (NaN ₃). Phosphate buffered saline
Form/Appearance	Recombinant Hepatitis B surface antigen AD
Concentration	0.5 mg/ml
Storage	+4°C DO NOT FREEZE
Note	For research use only
Clonality	Recombinant
Purity	Rec. Protein
Expiration Date	6 months from date of receipt.