

Product Datasheet

Bacterial CHLAMYDIA TRACHOMATIS protein (orb435897)

Description	Native bacterial CHLAMYDIA TRACHOMATIS protein
Reactivity	Bacteria
Tested Applications	ELISA, IA
Target	CHLAMYDIA TRACHOMATIS
Preservatives	None present. Phosphate buffered saline
Form/Appearance	Inactivated Chlamydia trachomatis - liquid
Concentration	4.9 mg/ml
Storage	-70°C (ship frozen)
Note	For research use only
Purity	Inactivated Pathogen
Expiration Date	6 months from date of receipt.