

## Product Datasheet

### Viral HERPES SIMPLEX VIRUS 1 protein (orb435901)

|                            |                                             |
|----------------------------|---------------------------------------------|
| <b>Description</b>         | Native viral HERPES SIMPLEX VIRUS 1 protein |
| <b>Reactivity</b>          | Virus                                       |
| <b>Tested Applications</b> | ELISA                                       |
| <b>Target</b>              | HERPES SIMPLEX VIRUS 1                      |
| <b>Preservatives</b>       | None present. Tissue culture medium         |
| <b>Form/Appearance</b>     | Inactivated Herpes simplex virus 1 - liquid |
| <b>Concentration</b>       | 2.14 mg/ml                                  |
| <b>Storage</b>             | -70°C (ship frozen)                         |
| <b>Note</b>                | For research use only                       |
| <b>Purity</b>              | Inactivated Pathogen                        |
| <b>Expiration Date</b>     | 6 months from date of receipt.              |