

Product Datasheet

Viral INFLUENZA B ANTIGEN protein (orb435903)

Description	Native viral INFLUENZA B ANTIGEN protein
Reactivity	Virus
Tested Applications	ELISA
Target	INFLUENZA B ANTIGEN
Preservatives	None present. Eagle's Minimum Essential Medium
Form/Appearance	Inactivated Influenza B antigen - liquid
Concentration	1.75 mg/ml
Storage	-70°C (ship frozen)
Note	For research use only
Purity	Inactivated Pathogen
Expiration Date	6 months from date of receipt.