

Product Datasheet

Viral MEASLES VIRUS protein (orb435904)

Description	Native viral MEASLES VIRUS protein
Reactivity	Virus
Tested Applications	ELISA
Target	MEASLES VIRUS
Preservatives	None present. Minimum Essential Medium
Form/Appearance	Inactivated Measles virus - liquid
Concentration	2.26 mg/ml.
Storage	-70°C (ship frozen)
Note	For research use only
Purity	Inactivated Pathogen
Expiration Date	6 months from date of receipt.