

Product Datasheet

Viral MUMPS VIRUS protein (orb435905)

Description	Native viral MUMPS VIRUS protein
Reactivity	Virus
Tested Applications	ELISA
Target	MUMPS VIRUS
Preservatives	None present. Culture medium 1969B
Form/Appearance	Inactivated Mumps virus - liquid
Concentration	2.35 mg/ml.
Storage	-70°C (ship frozen)
Note	For research use only
Purity	Inactivated Pathogen
Expiration Date	6 months from date of receipt.