


Product Datasheet

HLA-G antibody (FITC) (orb43591)

Description	Mouse Monoclonal to HLA-G.
Reactivity	Human
Conjugation	FITC
Tested Applications	FC
Immunogen	C-terminal amino acid sequence (22-mer) of soluble HLA-G5 and HLA-G6 proteins coupled to ovalbumin.
Target	HLA-G
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Excellent for intracellular staining; recommended dilution: 1-5 µg/ml.
Isotype	Mouse IgG1
Clonality	Monoclonal
Purity	Purified antibody is conjugated with fluorescein isothiocyanate (FITC) under optimum conditions and unconjugated antibody and free fluorochrome are removed by size-exclusion chromatography.
Uniprot ID	P17693
Entrez	3135
Dilution Range	Flow cytometry: Excellent for intracellular staining; recommended dilution: 1-5 µg/ml.
Expiration Date	12 months from date of receipt.



Flow
cytometry
analysis
(intracellular
S...