

Product Datasheet

Viral INFLUENZA A H1N1 ANTIGEN protein (orb435910)

Description	Native viral INFLUENZA A H1N1 ANTIGEN protein
Reactivity	Virus
Tested Applications	ELISA, Hemagglutination
Target	INFLUENZA A H1N1 ANTIGEN
Preservatives	None present. Phosphate buffered saline
Form/Appearance	Inactivated Influenza A H1N1 antigen - liquid
Concentration	2.2 mg/ml
Storage	-70°C (ship frozen)
Note	For research use only
Purity	Inactivated Pathogen
Expiration Date	6 months from date of receipt.