

## Product Datasheet

### Viral DENGUE VIRUS TYPE 1 NS1 ANTIGEN protein (orb435914)

|                            |                                                                     |
|----------------------------|---------------------------------------------------------------------|
| <b>Description</b>         | Recombinant of viral DENGUE VIRUS TYPE 1 NS1 ANTIGEN protein        |
| <b>Reactivity</b>          | Virus                                                               |
| <b>Tested Applications</b> | ELISA                                                               |
| <b>Target</b>              | DENGUE VIRUS TYPE 1 NS1 ANTIGEN                                     |
| <b>Preservatives</b>       | None present. Dulbecco's phosphate buffered saline                  |
| <b>Form/Appearance</b>     | Purified recombinant protein with a C-terminal 6 x His-tag - liquid |
| <b>Concentration</b>       | 0.46 mg/ml                                                          |
| <b>Storage</b>             | -70°C (ship frozen)                                                 |
| <b>Note</b>                | For research use only                                               |
| <b>Clonality</b>           | Recombinant                                                         |
| <b>Purity</b>              | Rec. Protein                                                        |
| <b>Expiration Date</b>     | 6 months from date of receipt.                                      |