

Product Datasheet

Viral DENGUE VIRUS TYPE 3 NS1 ANTIGEN protein (orb435918)

Description	Recombinant of viral DENGUE VIRUS TYPE 3 NS1 ANTIGEN protein
Reactivity	Virus
Tested Applications	ELISA
Target	DENGUE VIRUS TYPE 3 NS1 ANTIGEN
Preservatives	Purified recombinant protein with a C-terminal 6 x His-tag - liquid
Form/Appearance	Purified recombinant protein with a C-terminal 6 x His-tag - liquid
Concentration	0.5 mg/ml
Storage	-70°C (ship frozen)
Note	For research use only
Clonality	Recombinant
Purity	Rec. Protein
Expiration Date	6 months from date of receipt.