

Product Datasheet

Viral DENGUE VIRUS TYPE 4 NS1 ANTIGEN protein (orb435920)

Description	Recombinant of viral DENGUE VIRUS TYPE 4 NS1 ANTIGEN protein
Reactivity	Virus
Tested Applications	ELISA
Target	DENGUE VIRUS TYPE 4 NS1 ANTIGEN
Preservatives	None present. Tris buffered saline
Form/Appearance	Purified recombinant protein with a C-terminal His-tag - liquid
Concentration	0.54 mg/ml
Storage	-70°C (ship frozen)
Note	For research use only
Clonality	Recombinant
Purity	Rec. Protein
Expiration Date	6 months from date of receipt.