

Product Datasheet

Viral YELLOW FEVER VIRUS protein (orb435923)

Description	Recombinant of viral YELLOW FEVER VIRUS protein
Reactivity	Virus
Tested Applications	ELISA
Target	YELLOW FEVER VIRUS
Preservatives	Purified recombinant protein - liquid
Form/Appearance	Purified recombinant protein - liquid
Concentration	0.48 mg/ml
Storage	-70°C (ship frozen)
Note	For research use only
Clonality	Recombinant
Purity	Rec. Protein
Expiration Date	6 months from date of receipt.