

Product Datasheet

Mouse IL10 protein (orb435930)

Description	Recombinant of mouse IL10 protein
Reactivity	Mouse
Tested Applications	ELISA, FA, WB
Target	IL-10
Preservatives	Purified recombinant protein expressed in E. coli - lyophilized
Form/Appearance	Purified recombinant protein expressed in E. coli - lyophilized
Concentration	0.1mg/ml
Storage	-20°C only (ship +4°C)
Note	For research use only
Clonality	Recombinant
Purity	Rec. Protein
Uniprot ID	P18893
Expiration Date	6 months from date of receipt.