

Product Datasheet

Mouse MIP-1 ALPHA protein (orb435936)

Description	Recombinant of mouse MIP-1 ALPHA protein
Reactivity	Mouse
Tested Applications	ELISA, FA, WB
Target	MIP-1 ALPHA
Preservatives	None present
Form/Appearance	Purified recombinant protein - lyophilized
Concentration	0.1 mg/ml
Storage	-20°C after reconstitution
Note	For research use only
Clonality	Recombinant
Purity	Rec. Protein
Expiration Date	6 months from date of receipt.