

Product Datasheet

Mouse MIP-2 protein (orb435947)

Description	Recombinant of mouse MIP-2 protein
Reactivity	Mouse
Tested Applications	ELISA, FA, WB
Target	MIP-2
Preservatives	None present. 10mM Acetic Acid
Form/Appearance	Purified Recombinant Protein - lyophilized
Concentration	0.1mg/ml
Storage	-20°C after reconstitution
Note	For research use only
Clonality	Recombinant
Purity	Rec. Protein
Expiration Date	6 months from date of receipt.