

Product Datasheet

Pig MIP-1 BETA protein (orb435965)

Description	Recombinant of pig MIP-1 BETA protein
Reactivity	Porcine
Tested Applications	ELISA, WB
Target	MIP-1 BETA
Preservatives	10% Trehalose. Phosphate buffered saline
Form/Appearance	Purified recombinant protein - lyophilized
Storage	-20°C only (ship +4°C)
Note	For research use only
Clonality	Recombinant
Purity	Rec. Protein
Expiration Date	6 months from date of receipt.