

Product Datasheet

Rabbit IL1 ALPHA protein (orb435980)

Description	Recombinant of rabbit IL1 ALPHA protein
Reactivity	Rabbit
Tested Applications	ELISA, WB
Target	IL-1 ALPHA
Preservatives	10% Trehalose. Phosphate buffered saline
Form/Appearance	Purified recombinant protein - lyophilized
Storage	-20°C only (ship +4°C)
Note	For research use only
Clonality	Recombinant
Purity	Rec. Protein
Expiration Date	6 months from date of receipt.