

## Product Datasheet

### Rabbit IL8 protein (orb435983)

|                            |                                            |
|----------------------------|--------------------------------------------|
| <b>Description</b>         | Recombinant of rabbit IL8 protein          |
| <b>Reactivity</b>          | Rabbit                                     |
| <b>Tested Applications</b> | ELISA, WB                                  |
| <b>Target</b>              | IL-8                                       |
| <b>Preservatives</b>       | 10% Trehalose. Phosphate buffered saline   |
| <b>Form/Appearance</b>     | Purified recombinant protein - lyophilized |
| <b>Storage</b>             | -20°C only (ship +4°C)                     |
| <b>Note</b>                | For research use only                      |
| <b>Clonality</b>           | Recombinant                                |
| <b>Purity</b>              | Rec. Protein                               |
| <b>Expiration Date</b>     | 6 months from date of receipt.             |