

Product Datasheet

Rabbit EXODUS-2 protein (orb435986)

Description	Recombinant of rabbit EXODUS-2 protein
Reactivity	Rabbit
Tested Applications	ELISA, WB
Target	EXODUS-2
Preservatives	10% Trehalose. Phosphate buffered saline
Form/Appearance	Purified recombinant protein - lyophilized
Storage	-20°C only (ship +4°C)
Note	For research use only
Clonality	Recombinant
Purity	Rec. Protein
Expiration Date	6 months from date of receipt.