

Product Datasheet

Rat IL6 protein (orb435991)

Description	Recombinant of rat IL6 protein
Reactivity	Rat
Tested Applications	ELISA, FA, WB
Target	IL-6
Preservatives	None present. TRIS buffered saline0.1% SDS
Form/Appearance	Purified recombinant protein expressed in E. coli - lyophilized.
Concentration	0.1 mg/ml
Storage	-20°C after reconstitution
Note	For research use only
Clonality	Recombinant
Purity	Rec. Protein
Expiration Date	6 months from date of receipt.