

Product Datasheet

Rat IL10 protein (orb435996)

Description	Recombinant of rat IL10 protein
Reactivity	Rat
Tested Applications	ELISA, FA, WB
Target	IL-10
Preservatives	None present. Sodium Phosphate
Form/Appearance	Purified recombinant protein - lyophilized
Concentration	1.0 mg/ml
Storage	-20°C after reconstitution
Note	For research use only
Clonality	Recombinant
Purity	Rec. Protein
Expiration Date	6 months from date of receipt.