

Product Datasheet

STREPTAVIDIN protein (APC) (orb436001)

Description	STREPTAVIDIN protein (APC)
Conjugation	APC
Tested Applications	FC
Target	STREPTAVIDIN
Preservatives	0.09% Sodium Azide 1% Bovine Serum Albumin 5% Sucrose. Phosphate buffered saline
Form/Appearance	Streptavidin conjugated to Allophycocyanin - lyophilized
Concentration	0.1mg/ml
Storage	+4°C DO NOT FREEZE
Note	For research use only
Expiration Date	6 months from date of receipt.