

Product Datasheet

STREPTAVIDIN protein (FITC) (orb436002)

Description	STREPTAVIDIN protein (FITC)
Conjugation	FITC
Tested Applications	FC, IH-Fr, IHC-P
Target	STREPTAVIDIN
Preservatives	0.09% Sodium Azide (NaN ₃) 1% Bovine Serum Albumin. Phosphate buffered saline
Form/Appearance	Streptavidin conjugated to Fluorescein Isothiocyanate (Isomer 1) - liquid
Concentration	1 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Expiration Date	6 months from date of receipt.