

Product Datasheet

STREPTAVIDIN protein (RPE) (orb436003)

Description	STREPTAVIDIN protein (RPE)
Conjugation	RPE
Tested Applications	FC
Target	STREPTAVIDIN
Preservatives	0.09% Sodium Azide 1% Bovine Serum Albumin. Phosphate buffered saline
Form/Appearance	Streptavidin conjugated to R-Phycoerythrin - lyophilized
Storage	+4°C DO NOT FREEZE
Note	For research use only
Expiration Date	6 months from date of receipt.