

Product Datasheet

STREPTAVIDIN protein (AP) (orb436005)

Description	STREPTAVIDIN protein (AP)
Conjugation	AP
Tested Applications	ELISA, IH-Fr, IHC-P, WB
Target	STREPTAVIDIN
Preservatives	0.05% Sodium Azide (NaN ₃) 1% Bovine Serum Albumin. 0.05M TRIS Chloride 0.15M NaCl 0.001M MgCl ₂ 0.0001M ZnCl ₂ 50% (v/v) Glycerol; pH 8.0
Form/Appearance	Streptavidin conjugated to alkaline phosphatase - liquid
Concentration	1.0 mg/ml
Storage	+4°C DO NOT FREEZE
Note	For research use only
Expiration Date	6 months from date of receipt.