

Product Datasheet

IgG antibody (HRP) (orb436162)

Description	Goat polyclonal antibody against IgG conjugated to HRP
Species/Host	Goat
Reactivity	Human
Conjugation	HRP
Tested Applications	ELISA, WB
Immunogen	Human serum.
Target	IgG
Preservatives	0.01%Thiomersal 1%Bovine Serum Albumin. Phosphate buffered saline
Form/Appearance	Purified IgG conjugated to Horseradish Peroxidase (HRP) - liquid
Concentration	0.5 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Application notes	F(ab') ₂
Isotype	Polyclonal IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.