

Product Datasheet

IgM antibody (APC) (orb436165)

Description	Goat polyclonal antibody against IgM conjugated to APC
Species/Host	Goat
Reactivity	Mouse
Conjugation	APC
Tested Applications	ELISA, FACS, ICC, IF
Immunogen	Mouse IgM paraproteins.
Target	IgM
Preservatives	0.1% Sodium Azide (NaN ₃)Stabilizing agent (sucrose). Phosphate buffered saline
Form/Appearance	Purified Ig conjugated to Allophycocyanin (APC) - liquid
Concentration	0.5 mg/ml
Storage	+4°C DO NOT FREEZE
Note	For research use only
Isotype	Polyclonal IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.