

Product Datasheet

IgG antibody (AP) (orb436179)

Description	Rabbit polyclonal antibody against IgG conjugated to AP
Species/Host	Rabbit
Reactivity	Gallus
Conjugation	AP
Tested Applications	ELISA, IH-Fr, WB
Immunogen	Highly purified, whole molecule chicken IgG.
Target	IgG
Preservatives	0.1% Sodium Azide (NaN ₃)1% Bovine Serum Albumin. 0.05M Tris chloride 0.15NaCl 0.001M MgCl ₂ 0.0001M ZnCl ₂ 50% (v/v) Glycerol; pH 8.0
Form/Appearance	Purified IgG conjugated to Alkaline Phosphatase - liquid
Concentration	1.0mg/ml
Storage	+4°C DO NOT FREEZE
Note	For research use only
Application notes	(H/L)
Isotype	Polyclonal IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.